


REAL[®] 6100 - Technical data and measurements



DATA	STANDARD MEASUREMENTS
Top speed4,5 km/h	Total width650 mm
Turning space 180°930 mm	Length800 mm
Driving wheel.....Both center wheels	Wheel front/back125mm
Driving time per chargeapprox. 8 hours	Driving wheel.....210 mm
Obstacle ability40 mm	Seat height490-770, 440-640 mm
Weight74,5 kg incl. batt.	Leg support length370-530 mm
Charging timeapprox. 6-8 hours	Seat width330-480 mm
Charger tested and approved by	Seat depth350-480 mm
Mercado Medic ABMentzer 3A	Seat angle-2°
Batteries tested and approved by	Back support height ..360-420 mm
Mercado Medic ABSonnenschein 2pcs 12V 25 Ah	Back support angle-15-+45°

PROGRAMMING DEVICE

Programmable functions:

- Forward speed	- Turning speed
- Forward acceleration	- Turning acceleration
- Forward deceleration	- Turning deceleration
- Reverse speed	- Joystick sensitivity
- Reverse acceleration	- Reversed joystick function
- Reverse deceleration	

MERCADO MEDIC
 Mercado Medic AB
 Box 1074
 181 22 Lidingö
 Sweden

Tel: +46-(0)-8 555 143 00
 Fax: +46-(0)-8 555 143 99
 E-mail: info@mercado.se

Mercado Medic Export Office:
 Nobelweg 22
 3899 BN Zeewolde
 The Netherlands

Tel: +31-(0)-36 521 99 95
 Fax: +31-(0)-36 521 99 90
 E-mail: export@mercado.se

