


REAL[®] 6100 - Technical data and measurements



DATA

Top speed	4,5 km/h
Turning space 180°	930 mm
Driving wheel	Both center wheels
Driving time per charge	approx. 8 hours
Obstacle ability	40 mm
Weight	74,5 kg incl. batt.
Charging time	approx. 6-8 hours
Charger tested and approved by	
Mercado Medic AB	Mentzer 3A
Batteries tested and approved by	
Mercado Medic AB	Sonnenschein 2pcs 12V 25 Ah


STANDARD MEASUREMENTS

Total width	650 mm
Length	800 mm
Wheel front/back	125mm
Driving wheel	210 mm
Seat height	490-770, 440-640 mm
Leg support length	370-530 mm
Seat width	330-480 mm
Seat depth	350-480 mm
Seat angle	-2°
Back support height	360-420 mm
Back support angle	-15+45°

PROGRAMMING DEVICE

Programmable functions:

- | | |
|------------------------|------------------------------|
| - Forward speed | - Turning speed |
| - Forward acceleration | - Turning acceleration |
| - Forward deceleration | - Turning deceleration |
| - Reverse speed | - Joystick sensitivity |
| - Reverse acceleration | - Reversed joystick function |
| - Reverse deceleration | |



Mercado Medic AB
Box 1074
181 22 Lidingö
Sweden

Tel: +46-(0)-8 555 143 00
Fax: +46-(0)-8 555 143 99
E-mail: info@mercado.se

Mercado Medic Export Office:
Nobelweg 22
3899 BN Zeewolde
The Netherlands

Tel: +31-(0)-36 521 99 95
Fax: +31-(0)-36 521 99 90
E-mail: export@mercado.se

