

Try out/orderform REAL Mobil 6100 ERGOMEDIC

A	B	C	D	E	F	G	H	Upholstery
6	1							


 Offer request

 Order

Delivery address		Invoice address	
Tried out by		Tried out for/marketing	
Telephone	E-mail	Purchasers order number	Date

Basic equipment - Always fill in! AB Model 61 <input type="checkbox"/> REAL MOBIL 6100		Alternatives and accessories - fill in if you want to change or add something!	
C Seat angle adjustment 0 <input type="checkbox"/> Without adjustment 2 <input type="checkbox"/> Crank adjustment 3 <input type="checkbox"/> Gasspring adjustment at seat 4 <input type="checkbox"/> Gasspring adjustment at armrest 5 <input type="checkbox"/> Electric adjustment 6 <input type="checkbox"/> Prepared for seat angle adjustment		Programming unit <input type="checkbox"/> Yes <input type="checkbox"/> No Mark Yes if programming unit is needed	
D Low back rest 0 <input type="checkbox"/> Standard back rest 1 <input type="checkbox"/> Medic back rest 2 <input type="checkbox"/> Comfort back rest* 3 <input type="checkbox"/> Electric back rest*		Upholstery Plush: <input type="checkbox"/> 41 bordeaux <input type="checkbox"/> 42 brown <input type="checkbox"/> 43 green <input type="checkbox"/> 44 marine <input type="checkbox"/> 45 REAL-blue <input type="checkbox"/> 46 grey <input type="checkbox"/> 48 black Patterned fabric Adult: <input type="checkbox"/> 71 red <input type="checkbox"/> 75 blue <input type="checkbox"/> 78 black Patterned fabric Child: <input type="checkbox"/> 81 red <input type="checkbox"/> 85 blue <input type="checkbox"/> 89 yellow Skai: <input type="checkbox"/> 31 red <input type="checkbox"/> 32 brown <input type="checkbox"/> 35 blue <input type="checkbox"/> 36 grey <input type="checkbox"/> 38 black	
D High back rest 5 <input type="checkbox"/> Standard back rest 6 <input type="checkbox"/> Medic back rest <input type="checkbox"/> Double Medic 7 <input type="checkbox"/> Comfort back rest* 8 <input type="checkbox"/> Electric back rest* <i>* Only for seat 43 cm depth or more</i>		Incontinence upholstery (add cover) <input type="checkbox"/> Seat <input type="checkbox"/> Backrest <input type="checkbox"/> Leg divider Covers: <input type="checkbox"/> Seat <input type="checkbox"/> Backrest <input type="checkbox"/> Others according to below	
E Leg support 0 Standard with foldable foot rest <input type="checkbox"/> 30-38 cm <input type="checkbox"/> 37-53 cm 9 <input type="checkbox"/> Alternative (choose to the right)		Alternative cotton covers: <input type="checkbox"/> Child/yellow <input type="checkbox"/> Child/blue <input type="checkbox"/> Incontinence cover, grey	
F Seat height seat hight (angleadj.) 0 <input type="checkbox"/> Standard 49-77 cm 50-78 cm 1 <input type="checkbox"/> Low 44-64 cm 48-68 cm		Arm rests (foldable) <input type="checkbox"/> Black skai <input type="checkbox"/> Non foldable <input type="checkbox"/> Plush covered <input type="checkbox"/> Natural leather Length: <input type="checkbox"/> 35 cm <input type="checkbox"/> 30 cm <input type="checkbox"/> 25 cm	
G Seat "Slow Recovery" widthxdepth 0 <input type="checkbox"/> Adult 44x48 cm 1 <input type="checkbox"/> Adult 48x48 cm 2 <input type="checkbox"/> Adult 44x43 cm 3 <input type="checkbox"/> Adult/child 40x43 cm 4 <input type="checkbox"/> Adult/child 40x39 cm 5 <input type="checkbox"/> Child (7-12 years) 36x39 cm 6 <input type="checkbox"/> Child (3-7 years) 33x35 cm 7 <input type="checkbox"/> Child 36x43 cm 8 <input type="checkbox"/> Adult 40x48 cm		Backrest width x height in cm Armrest Low 38x26 extra high 42x43 35 cm Low 38x26 extra high 42x43 35 cm Low 38x26 extra high 42x43 35 cm Low 33x26 high 33x36 35 cm Low 33x26 high 33x36 30 cm Low 30x26 high 30x34 30 cm Low 27x22 high 27x32 30 cm Low 30x26 high 30x34 30 cm Low 33x26 high 33x36 35 cm	
H User's weight 0 <input type="checkbox"/> 55-90 kg 1 <input type="checkbox"/> < 55 kg User's weight >130 kg 2 <input type="checkbox"/> 90-130 kg contact Mercado		EXTRA	
Battery <input type="checkbox"/> Yes <input type="checkbox"/> No Battery charger <input type="checkbox"/> Yes <input type="checkbox"/> No		No. of chairs:	

Seat and back rest in ErgoMedic PLUS <input type="checkbox"/> Seat <input type="checkbox"/> Seat Slow Recovery <input type="checkbox"/> Back rest 40x20 cm <input type="checkbox"/> Back rest 30x34 cm <input type="checkbox"/> Back rest 35x43 cm <input type="checkbox"/> Back rest 42x43 cm <input type="checkbox"/> Back rest 46x43 cm <input type="checkbox"/> Back rest 35x43 cm with wedge 7 cm <input type="checkbox"/> Back rest 42x43 cm with wedge 7 cm <input type="checkbox"/> Back rest 46x43 cm with wedge 7 cm	
--	--

Placing of device Joystick <input type="checkbox"/> Right <input type="checkbox"/> Left

Amputation rests <input type="checkbox"/> Right <input type="checkbox"/> Left

Leg rest Medic with foot plate Leg rest undivided <input type="checkbox"/> S <input type="checkbox"/> M <input type="checkbox"/> L <input type="checkbox"/> XL Leg rest divided <input type="checkbox"/> S <input type="checkbox"/> M <input type="checkbox"/> L <input type="checkbox"/> XL

Electric leg rest <input type="checkbox"/> Right <input type="checkbox"/> Left <input type="checkbox"/> Right + Left
--

Adapter Cross for leg support <input type="checkbox"/> Right <input type="checkbox"/> Left <input type="checkbox"/> Right + Left
--

Trunk support (mounted on high backrest) <input type="checkbox"/> Standard <input type="checkbox"/> Rounded <input type="checkbox"/> Hight and depth adjustable
--

Hip belt (with attachment under the seat) <input type="checkbox"/> S, 51-71 <input type="checkbox"/> M, 63-83 <input type="checkbox"/> L, 83-103 cm

Knee support (for leg support Medic) <input type="checkbox"/> 15x8 cm black, pair

Calf support (for leg support Mobil/Medic) <input type="checkbox"/> 15x13 cm
--

Neck rests (mounted on high back) <input type="checkbox"/> Flat <input type="checkbox"/> Curved <input type="checkbox"/> ErgoMedic Plus

Side support adjustable in width (wxh)* <input type="checkbox"/> 28x8 <input type="checkbox"/> 24x14 <input type="checkbox"/> 30x16 <input type="checkbox"/> 15x8 <i>*Combine with non foldable armrests</i>

REAL 6100/0806